

Integrated NSK Motor Selection Guide

ASI Dental can integrate the various models of the NSK electric motor line up. The following provides a brief description of the various models and guidance on consideration for integration.



NLX Nano Electric Motor - NSK PN Y1001361

Highly Recommended for Cart Integration - This motor control features a compact panel design. It integrated internally flush with ASI panels. Two motors can be installed with two controls that can be mounted together in one panel.

NLZ Electric Motor - NSK PN Y1003735

Acceptable for Integration - This control panel can be integrated internally flush with ASI panels. The control and display do require a full ASI panel. So if you want two motors, will require a second panel. *The NSK NLZ is not recommended for side mounting in Classic Series Carts. The cart's compact footprint requires the motor's control panel to be oriented sideways rather than upright.*



NLZE Electric Motor with Endo Function - NSK PN Y1003736

Recommend for Reciprocating function if desired. Acceptable for Integration - This control panel can be integrated internally flush with ASI panels. The control and display do require a full ASI panel. So, if you want two motors, it will require a second panel. The reciprocating function can complement using a standard rotary motor control such as DTC without having to use a ProMark Model. ProMark motors do require a dual panel set up.

NLZ Pro Electric Motor – NSK PN Y1500670

NOT Recommended for Integration – This motor can run one to two motors. However, the control cannot be internally integrated and must be mounted on top of the ASI control panel or use a separate bracket on the side. The cord from the NLZ Pro will also be run outside from it and then back inside the ASI cart.



ASIdental.com | 844.880.3636

8811 American Way, Suite 130 | Englewood, CO 80112
Advanced Endodontic Systems® | Advanced Dental Systems®

 **ASI**
Dental Specialties

08-05-26revB